

# NEW – NEW – NEW THE BITE SPLINT WITH THERMO-EFFECT



excited  
patients

pleased  
dentists

relaxed  
technicians

**Nothing presses, nothing straints**

**the milled PREMIOtemp CLEAR FLEX bite splint is just super comfortable to wear.**

## PREMIOtemp CLEAR FLEX

the **NEW** standard for milled **bite splints!**

- **PREMIOtemp CLEAR FLEX is the new** and very innovative milling blank for **CAD/CAM milled nightguards** to our primotec digital products portfolio.
- The CLEAR FLEX material is a **transparent PMMA acrylic with thermo-effect**. This effect procures that the milled bite splint, when warmed up, automatically adjusts to the dental arch in the patients mouth, without losing the materials surface hardness.
- The thermo-effect is achieved, when the milled PREMIOtemp CLEAR FLEX nightguard is **pre-heated in 40°C to 50°C warm water** prior to being inserted. This way, during insertion, the bite splint perfectly adapts itself in the patient's mouth.

>>>>

# CLEAR FLEX



# THE NEW STANDARD FOR MILLED BITE SPLINTS



0  
1  
2  
3  
4  
5  
6  
7  
8  
9

The result is a milled nightguard that **perfectly fits the dental arch without any tension**. Due to the surface hardness of the material, occlusal contacts and canine- as well as anterior guidance remain where they had been set in the lab.

**Especially comfortable for the dentist and the lab** is the fact, that PREMIOtemp CLEAR FLEX, due to the described thermo-effect, **can compensate for moderate inaccuracies of the impression**. This way, the material is especially appropriate **for intra oral scanning of the complete dental arch**.

## CONCLUSION:

- **excited patients**
- **pleased dentists**
- **relaxed technicians**

**due to milled PREMIOtemp CLEAR FLEX nightguards from primotec**

**Contact: +49(0) 61 72 - 99 770 - 0**

**or: [primotec@primogroup.de](mailto:primotec@primogroup.de)**